



SP V 3DH12E RO-D



SP V 3DC12B RO-D



SP V 3DC12V RO-DV

Product Features

- oil free
- maintenance free
- reduced vibrations
- low noise
- small size
- lightweight
- minimal current consumption


Applications

- Blood Pressure Monitoring




Designation	SP V 3DH12E RO-D	SP V 3DC12B RO-D	SP V 3DC12V RO-DV
Voltage	3 Vdc	3 Vdc	3 Vdc
optional Voltage	-	-	6 Vdc
Art.-No.	7k73018	7k73016	7k73029
Pneumatic Performance			
Free Flow [l/min]	0,70	0,66	0,38
flow at:			
100 mbar [l/min]	0,51	0,54	0,32
max. Pressure [mbar]	620	830	680
flow at:			
-100 mbar [l/min]	-	-	0,28
max. Vacuum [mbar]	-	-	-400
Motor / Power Consumption			
Motor Type		Iron Core	
max. Current [mA]	250	360	350
Construction			
Pump Body		ABS / POM	
Diaphragm		EPDM	
Valves		NBR	
Weight [g]	11	23	23

Performance Curves






SP V 3DH12E RO-D



SP V 3DC12B RO-D



SP V 3DC12V RO-DV

